



## CERTIFICATE OF ANALYSIS

This "Certificate of Analysis" represents a precleaned product that has been prepared in accordance with Performance-Based specifications. This product meets or exceeds analyte specifications established in the U. S. EPA OSWER Directive 9240.0-05A "Specification and Guidance for Contaminant-free Sample Containers" for use in Superfund and other Hazardous waste programs.

### Group 3 Volatile Organic Compounds (VOCs)

| Analyte                     | RL ug/L | Analyte                    | RL ug/L | Analyte                   | RL ug/L    |
|-----------------------------|---------|----------------------------|---------|---------------------------|------------|
| Acetone                     | 5.0 U   | trans-1,2-Dichloroethylene | 0.50 U  | n-Propylbenzene           | 0.50 U     |
| Benzene                     | 0.50 U  | 1,2-Dichloropropane        | 0.50 U  | Styrene                   | 0.50 U     |
| Bromobenzene                | 0.50 U  | 1,3-Dichloropropane        | 0.50 U  | Tert-Amyl Alcohol         | 50 U       |
| Bromochloromethane          | 0.50 U  | 2,2-Dichloropropane        | 0.50 U  | Tert-Amyl Methyl Ether    | 1.0 U      |
| Bromodichloromethane        | 0.50 U  | 1,1-Dichloropropene        | 0.50 U  | Tert-Butyl Alcohol        | 50 U       |
| Bromoform                   | 0.50 U  | cis-1,3-Dichloropropene    | 0.50 U  | 1,1,1,2-Tetrachloroethane | 0.50 U     |
| n-Butylbenzene              | 0.50 U  | trans-1,3-Dichloropropene  | 0.50 U  | 1,1,2,2-Tetrachloroethane | 0.50 U     |
| sec-Butylbenzene            | 0.50 U  | Di-Isopropyl Ether         | 0.50 U  | Tetrachloroethylene       | 0.50 U     |
| tert-Butylbenzene           | 0.50 U  | Ethyl Acetate              | 0.50 U  | Toluene                   | 0.50 U     |
| Carbon Disulfide            | 0.50 U  | Ethyl Alcohol              | 50 U    | 1,2,3-Trichlorobenzene    | 0.50 U     |
| Carbon Tetrachloride        | 0.50 U  | Ethylbenzene               | 0.50 U  | 1,2,4-Trichlorobenzene    | 0.50 U     |
| Chlorobenzene               | 0.50 U  | Ethyl Tert-Butyl Ether     | 0.50 U  | 1,1,1-Trichloroethane     | 0.50 U     |
| Chloroethane                | 0.50 U  | 2-Hexanone                 | 1.0 U   | 1,1,2-Trichloroethane     | 0.50 U     |
| Chloroform                  | 0.50 U  | Hexachlorobutadiene        | 0.50 U  | Trichloroethylene         | 0.50 U     |
| 2-Chlorotoluene (ortho)     | 0.50 U  | Isopropyl Alcohol          | 5.0 U   | Trichlorofluoromethane    | 0.50 U     |
| 4-Chlorotoluene (para)      | 0.50 U  | Isopropylbenzene           | 0.50 U  | 1,2,3-Trichloropropane    | 0.50 U     |
| Dibromochloromethane        | 0.50 U  | p-Isopropyltoluene         | 0.50 U  | 1,2,4-Trimethylbenzene    | 0.50 U     |
| 1,2-Dibromo-3-chloropropane | 0.50 U  | Methyl Acetate             | 0.50 U  | 1,3,5-Trimethylbenzene    | 0.50 U     |
| 1,2-Dibromoethane           | 0.50 U  | Methyl Bromide             | 0.50 U  | Vinyl Acetate             | 0.50 U     |
| Dichlorodifluoromethane     | 0.50 U  | Methyl Chloride            | 0.50 U  | Vinyl Chloride            | 0.50 U     |
| 1,3-Dichlorobenzene (meta)  | 0.50 U  | Methylene Bromide          | 0.50 U  | m,p-Xylene                | 1.0 U      |
| 1,2-Dichlorobenzene (ortho) | 0.50 U  | Methylene Chloride         | 0.50 U  | o-Xylene                  | 0.50 U     |
| 1,4-Dichlorobenzene (para)  | 0.50 U  | Methyl Ethyl Ketone        | 5.0 U   | Xylene (total)            | 0.50 U     |
| 1,1-Dichloroethane          | 0.50 U  | 4-Methyl-2-pentanone       | 0.50 U  | GRO (8015 C6-C10)         | 100 U      |
| 1,2-Dichloroethane          | 0.50 U  | Methyl Tert Butyl Ether    | 0.50 U  | TOC                       | 1.0 U mg/L |
| 1,1-Dichloroethylene        | 0.50 U  | Naphthalene                | 0.50 U  |                           |            |
| cis-1,2-Dichloroethylene    | 0.50 U  |                            |         |                           |            |

#### NOTES:

- RL = Reporting Limit.
- U = The analyte was analyzed for but not detected above the Reporting Limit.
- Bottles are Type III Soda Lime and vials are Type I Borosilicate.
- Vial Storage: Store at 85F or 29.4C. Keep away from organic vapors.
- Solid-top caps feature fluoropolymer resin liners. Open-top caps feature ultrasonically bonded 3.1mm (1/8") fluoropolymer resin/silicone septa.

This "Certificate of Analysis" is provided for your records and is used to facilitate any required correspondences as needed.

This case is preserved. These containers are barcoded.

Each container contains: 25 mg Ascorbic Acid

#### Level: Quality Assured (QA) (X)

Glassware / Plasticware received full Quality Assurance and Quality Control treatment. Containers, liners, and closures as applicable, are cleaned according to EPA recommended procedures and validated through a third party (NELAP) testing Laboratory. Each case of containers is custody sealed and labeled for traceability by Lot Number.

Part Number: QLAV89080066

Lot Number: 072026-3BNB

VWR Part No.: 89094-180

Date Product Prepared: 7/23/2026

Item Description: 40mL Amber Borosilicate Vial

Chemical Lot Number: A0467131

Protocol: B

Level: QA

Chemical Expiry Date: 11/04/29

Group: 3 (applies)

Chief Executive Officer  
Manufactured for VWR International

B

40mL Amber Borosilicate Vial



## CERTIFICATE OF ANALYSIS

### Ascorbic Acid

This "Certificate of Analysis" represents a chemical preservative that has been prepared in accordance with Performance-Based specifications. This product meets or exceeds analyte specifications established in the U. S. EPA OSWER Directive 9240.0-05A "Specification and Guidance for Contaminant-free Sample Containers" for use in Superfund and other Hazardous waste programs.

The Ascorbic Acid used in preservation procedures is tested by a 3rd Party NELAC Laboratory prior to use for determination of potential volatile contaminants based on the list below.

### Group 3 Volatile Organic Compounds (VOCs)

| <u>Analyte</u>              | <u>RL ug/L</u> | <u>Analyte</u>             | <u>RL ug/L</u> | <u>Analyte</u>            | <u>RL ug/L</u> |
|-----------------------------|----------------|----------------------------|----------------|---------------------------|----------------|
| Acetone                     | 5.0 U          | trans-1,2-Dichloroethylene | 0.50 U         | n-Propylbenzene           | 0.50 U         |
| Benzene                     | 0.50 U         | 1,2-Dichloropropane        | 0.50 U         | Styrene                   | 0.50 U         |
| Bromobenzene                | 0.50 U         | 1,3-Dichloropropane        | 0.50 U         | Tert-Amyl Alcohol         | 50 U           |
| Bromochloromethane          | 0.50 U         | 2,2-Dichloropropane        | 0.50 U         | Tert-Amyl Methyl Ether    | 1.0 U          |
| Bromodichloromethane        | 0.50 U         | 1,1-Dichloropropene        | 0.50 U         | Tert-Butyl Alcohol        | 50 U           |
| Bromoform                   | 0.50 U         | cis-1,3-Dichloropropene    | 0.50 U         | 1,1,1,2-Tetrachloroethane | 0.50 U         |
| n-Butylbenzene              | 0.50 U         | trans-1,3-Dichloropropene  | 0.50 U         | 1,1,2,2-Tetrachloroethane | 0.50 U         |
| sec-Butylbenzene            | 0.50 U         | Di-Isopropyl Ether         | 0.50 U         | Tetrachloroethylene       | 0.50 U         |
| tert-Butylbenzene           | 0.50 U         | Ethyl Acetate              | 0.50 U         | Toluene                   | 0.50 U         |
| Carbon Disulfide            | 0.50 U         | Ethyl Alcohol              | 50 U           | 1,2,3-Trichlorobenzene    | 0.50 U         |
| Carbon Tetrachloride        | 0.50 U         | Ethylbenzene               | 0.50 U         | 1,2,4-Trichlorobenzene    | 0.50 U         |
| Chlorobenzene               | 0.50 U         | Ethyl Tert-Butyl Ether     | 0.50 U         | 1,1,1-Trichloroethane     | 0.50 U         |
| Chloroethane                | 0.50 U         | 2-Hexanone                 | 1.0 U          | 1,1,2-Trichloroethane     | 0.50 U         |
| Chloroform                  | 0.50 U         | Hexachlorobutadiene        | 0.50 U         | Trichloroethylene         | 0.50 U         |
| 2-Chlorotoluene (ortho)     | 0.50 U         | Isopropyl Alcohol          | 5.0 U          | Trichlorofluoromethane    | 0.50 U         |
| 4-Chlorotoluene (para)      | 0.50 U         | Isopropylbenzene           | 0.50 U         | 1,2,3-Trichloropropane    | 0.50 U         |
| Dibromochloromethane        | 0.50 U         | p-Isopropyltoluene         | 0.50 U         | 1,2,4-Trimethylbenzene    | 0.50 U         |
| 1,2-Dibromo-3-chloropropane | 0.50 U         | Methyl Acetate             | 0.50 U         | 1,3,5-Trimethylbenzene    | 0.50 U         |
| 1,2-Dibromoethane           | 0.50 U         | Methyl Bromide             | 0.50 U         | Vinyl Acetate             | 0.50 U         |
| Dichlorodifluoromethane     | 0.50 U         | Methyl Chloride            | 0.50 U         | Vinyl Chloride            | 0.50 U         |
| 1,3-Dichlorobenzene (meta)  | 0.50 U         | Methylene Bromide          | 0.50 U         | m,p-Xylene                | 1.0 U          |
| 1,2-Dichlorobenzene (ortho) | 0.50 U         | Methylene Chloride         | 0.50 U         | o-Xylene                  | 0.50 U         |
| 1,4-Dichlorobenzene (para)  | 0.50 U         | Methyl Ethyl Ketone        | 25 U           | Xylene (total)            | 0.50 U         |
| 1,1-Dichloroethane          | 0.50 U         | 4-Methyl-2-pentanone       | 0.50 U         | TICs                      | U*             |
| 1,2-Dichloroethane          | 0.50 U         | Methyl Tert Butyl Ether    | 0.50 U         | Sulfide                   | 0.05 U mg/L    |
| 1,1-Dichloroethylene        | 0.50 U         | Naphthalene                | 0.50 U         | Cyanide                   | 0.02 U mg/L    |
| cis-1,2-Dichloroethylene    | 0.50 U         |                            |                |                           |                |

#### NOTES:

- Reporting Limit (RL) = The lowest concentration standard analyzed which can be verified.
- U = The analyte was analyzed for but not detected above the Reporting Limit.
- U\* = No analytes were detected; No Reporting Limits for these analytes.
- Solid-top caps feature fluoropolymer resin liners. Open-top caps feature ultrasonically bonded 3.1mm (1/8") fluoropolymer resin/silicone septa.

The chemical added to this product may become altered depending upon storage conditions and elapsed time from the prepared date.

This "Certificate of Analysis" is provided for your records and is used to facilitate any required correspondences as needed.

Item Description: Ascorbic Acid (ACS Grade)

Lot Number: A0467131

Date Product Prepared: 07/23/26

Chemical Expiry Date: 11/04/29

Chief Executive Officer  
Manufactured for VWR International